FLXIFIT™ 407 NASAL MASK PATIENT INSTRUCTIONS FOR USE

THE FLXIFIT™ 407 NASAL MASK IS INTENDED FOR MULTIPLE PATIENT OR SINGLE PATIENT ADULT USE BY INDIVIDUALS WHO HAVE BEEN DIAGNOSED BY A PHYSICIAN AS REQUIRING CPAP OR BI-LEVEL VENTILATOR TREATMENT.

FITTING YOUR MASK

- Hold mask in hand, spread the Stretchgear™ Headgear (H) with your free hand, just enough to fit comfortably over your head.
- Lightly hold the mask against your face covering the nose, then slowly pull the Headgear over your head.
- Attach the Headgear to the Glider™ Strap (E).
- Pull the lower straps to fit snugly over the back of your neck. Gently tighten the straps to remove slack. Connect to CPAP/Bi-Level machine set as prescribed.
- If necessary, tighten the lower straps to eliminate any leaks (avoid over-tightening the straps).
- When you need to remove your mask, unhook the Stretchgear™ Headgear from the Glider™ Strap.

FITTING TIPS

- Forehead pad should rest gently on the forehead without top straps being pulled too tight.
- If there are any leaks after gently tightening the top and lower straps, pull the mask slightly out from the face while straps are still attached. This will allow Silicone Seal to inflate with air. Put the mask back onto the face.

CLEANING YOUR MASK

To ensure continuing safety and comfort for your FlexFit™ 407 Nasal Mask:

BEFORE FIRST USE:

- Hand-wash the Stretchgear™ Headgear (H) in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes, then rinse in fresh water and allow to dry out of direct sunlight.

BEFORE EACH USE:

- Inspect the mask for deterioration. Do not use if damaged.

AFTER EACH USE:

- Disassemble the Silicone Seal (C), Foam Cushion (D) and the Headgear (H) from the Mask Base (A).
- Wash the Silicone Seal and Mask Base in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes, then rinse in fresh water.
- Wipe the foam cushion with a damp cloth. (Do not immerse the Foam Cushion in water.)
- Leave all parts to dry out of direct sunlight before reassembling.

EVERY 7 DAYS:

- Hand-wash the Stretchgear™ Headgear (H) in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes, then rinse in fresh water and allow to dry out of direct sunlight.

CAUTIONS

- Do not leave the mask for longer than 10 minutes.
- Do not wash the mask in the dishwasher.
- Do not clean the mask with products containing alcohol, anti-bacterial agents, antiseptic, bleach, chlorine or moisturizer.
- Do not store the mask in direct sunlight.

The above actions may deteriorate or damage the mask and shorten its life. If your FlexFit™ 407 Nasal Mask weakens or cracks, discontinue use and seek replacement immediately.

ASSEMBLING YOUR MASK

Your FlexFit™ 407 Nasal Mask is fully assembled and ready to use. However, when you need to clean your mask follow these instructions for reassembly.

- Place Silicone Seal (C) face down onto flat surface. With thick side of Foam Cushion (D) facing down, place into Silicone Seal.
- Starting with the slot at the apex of the Silicone Seal, assemble onto Mask Base (A). Push firmly together.
- Lay the Stretchgear™ Headgear (H) flat with the loop side down. Position the assembled mask over the Headgear and attach the four straps into the corresponding slots in the Forehead Rest (F) and the Glider™ Strap (E). This can be done without undoing the VELCRO™ brand tabs by sliding the Headgear into the slots.

*VELCRO™ is a registered trademark of Velcro Industries BV.

WARRANTY AND DISPOSAL

WARRANTY STATEMENT

Fisher & Paykel Healthcare warrants that the mask (excluding foams and material diffusers), when used in accordance with its instructions for use, will be free from defects in workmanship and materials and will perform in accordance with Fisher & Paykel Healthcare's official published product specifications for a period of 90 days from the date of purchase by the user. This warranty is subject to the limitations and exceptions set out in detail in <http://www.fphcare.com/osa/productwarranty.asp>

DISPOSAL

This device does not contain any hazardous substances and may be disposed of with general waste.

WARNINGS

- This mask should be used only with CPAP or Bi-Level equipment recommended by your physician or respiratory therapist. A mask should not be used unless the CPAP or Bi-Level machine is turned on and operating properly.
- Explanation of Warning:** CPAP and Bi-Level machines are intended to be used with special masks which have exhalation ports to allow continuous flow of air out of the mask. When the CPAP/Bi-Level machine is turned on and functioning properly, new air from the CPAP/Bi-Level machine flushes the exhaled air out through the mask exhalation ports. However, when the CPAP/Bi-Level machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to asphyxiation.
- At low CPAP pressures the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Do not block the gas exhaust ports. Before use, always ensure that gas is flowing out the exhaust ports.
- Do not use the mask if vomiting or feeling nauseous.
- Discontinue using the mask if there is an allergic reaction to any part of the mask. Consult your physician if this occurs.
- If re-using the FlexFit™ 407 between different patients, refer to the Multi-Patient Cleaning Instructions for the FlexFit™ 407 Nasal Mask.

CAUTIONS

- US Federal Law restricts this device to sale by, or on order of, a physician.

The FlexFit™ 407 Nasal Mask and accessories are latex free. Many people find CPAP/Bi-Level therapy causes dryness of the nose and throat. Heated humidification can substantially reduce these symptoms. Contact your supplier and ask about ThermoSmart™ Technology.

FOR ADDITIONAL INFORMATION CONTACT YOUR LOCAL HOMECARE PROVIDER

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OXYGEN/PRESSURE PORT CONNECTOR

If pressure readings and/or additional oxygen are required, an Oxygen/Pressure Port Connector is available (Code 900HC452).

- Connect the large end of the Connector (T) to the Swivel (G) on the mask and the small end of the Connector to the main breathing tube (R).

- Open the Part Cap(s) (P) and place the Pressure and/or Oxygen Tube (S) firmly over the Part(s) (Q).

(The two parts are identical; either part may be used for oxygen or pressure measurement. Only remove the Part Cap(s) from parts that are being used.)

- Note:** At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen will vary depending on the pressure settings, patient's breathing pattern, mask selection, and the leak rate.


- Warning:** If oxygen is used with CPAP/Bi-Level machine, the oxygen flow must be turned off when the CPAP/Bi-Level machine is not operating.

Explanation of Warning: When the CPAP/Bi-Level device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP/Bi-Level machine enclosure. Oxygen accumulated in the CPAP/Bi-Level machine enclosure will create a risk of fire.

- Warning:** Do not smoke while using or in the vicinity of someone using oxygen with this mask.

- Warning:** Do not use Polyvinyl Chloride (PVC) pressure and/or oxygen tubing. Use silicone tubing instead.

Explanation of Warning: Contact with PVC can cause premature cracking of parts.


TECHNICAL SPECIFICATIONS

EXHAUST FLOW

Nasal masks for CPAP and Bi-Level machines have exhaust flow holes to expel the air that you exhale from the mask.

PRESSURE (cmH ₂ O)	3	5	7	9	11	13	15	17	19	21	23	25
FLOW (L/MIN)	21	24	28	31	36	39	41	44	47	50	52	54

RESISTANCE TO FLOW

Pressure drop through mask @ 50 l/min 0.1 ± 0.1 cmH₂O

Pressure drop through mask @ 100 l/min 0.4 ± 0.1 cmH₂O

MASK DEAD SPACE

65cc

OPERATING CONDITIONS

- Operating pressure is 3-25 cmH₂O
- Operating temperature range is 5-50°C (40-120°F)

SOUND

The Sound Power Level of the mask is 31 dBA, with uncertainty ±2.5 dBA. The Sound Pressure Level of the mask is 25 dBA, with uncertainty ±2.5 dBA.

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